

# Performance of ThyroSeq v3 Genomic Classifier in Fixed Cytology Smears



L. TOLINO<sup>1</sup>, M. NIKIFOROVA<sup>1</sup>, G. TRONCONE<sup>2</sup>, A. WALD<sup>1</sup>, Z. BALOCH<sup>3</sup>, Y. NIKIFOROV<sup>1</sup>

## CONTACT INFORMATION

Marina N. Nikiforova, MD, Department of Pathology, University of Pittsburgh, 3477 Euler Way, Rm 8031, Pittsburgh, PA 15213 (nikiforovamn@upmc.edu)

<sup>1</sup>Department of Pathology, University of Pittsburgh, Pittsburgh, Pennsylvania


<sup>2</sup>Department of Public Health, University of Naples Federico II, Naples, Italy

<sup>3</sup>Department of Pathology and Laboratory Medicine, Perelman School of Medicine, University of Pennsylvania, Philadelphia

## INTRODUCTION

ThyroSeq v3 Genomic Classifier (GC) is a molecular test used to improve diagnosis in thyroid nodules with indeterminate fine-needle aspiration (FNA) cytology and inform patient management.

### Scheme of ThyroSeq Test Flow



FNA samples collected into a dedicated preservative solution vial is a standard sample type for ThyroSeq testing. However, in some cases, routine FNA smears are the only specimen type available.

## AIMS


- To evaluate ThyroSeq v3 GC analytical performance in:
1. FNA smears prepared for cytopathologic interpretation
  2. Reference cytology slides with gene mutations at varying allelic frequencies

## METHODS

- ThyroSeq v3 GC analysis was performed on:
- 33 clinical routinely prepared Diff-Quik<sup>®</sup> or Papanicolaou (PAP) stained slides, with ~200 to 1000 cells per slide, from 14 nodules
  - Four reference cytology slides (Horizon Diagnostics) with various mutation frequency
- Results were compared to FNAs collected into preservative solution (ThyroSeqPreserve)

## RESULTS

### ThyroSeq Performance in Cytology Smears (n=33, 14 nodules)




- In cytology smears, adequate DNA results for mutations and copy number alterations were obtained in 93% (13/14) of samples and adequate RNA results for gene fusions and gene expression in 79% (11/14) of samples.
- Two smears with <200-300 cells showed failure of RNA analysis.
- One old cellular smear (6 yrs) failed RNA due degraded nucleic acids.

### Comparison of Test Informative Rate Using Cytology Smears vs. Dedicated Collection Tube

FNA Preparation	Informative Result Rate	Required Cellularity	Turnaround Time
Fixed cytology smears (Diff-Quik, Pap)	79-93%	200-300 cells	8-9 days
FNA in ThyroSeqPreserve	96%	200-300 cells	5-7 days

- The rate of informative results was higher in FNAs collected into preservative solution as compared to cytology smears (79-93% vs. 96%).
- The turnaround time was 1-3 days longer due to additional preparatory steps for cytology smears as compared to FNA samples collected into preservative solution.

### Accuracy of Detection of Mutations in Reference Cytology Smears from World-Wide Ring Trial Study<sup>1</sup>



- Unstained reference cytology slides (Horizon Diagnostics) used in a Worldwide Ring Trial Study<sup>1</sup> with *KRAS*, *NRAS*, *BRAF*, and *PIK3CA* mutations at various frequency (10%, 5%, 1%, and 0%) were used to evaluate the accuracy of mutation detection in cytology smears.
- ThyroSeq detected all mutations down to 5% and *BRAF* and *PIK3CA* mutations down to 1% allele frequency. No mutations were detected in the negative control slide.

## CONCLUSIONS

- The results of this study provide evidence that ThyroSeq v3 GC testing can be effectively performed using routinely prepared cytology smears
- Overall, the test informative results were obtained in 79-93% of cytology smears
- The minimum required cellularity of cytology smears is 200-300 cells
- Both PAP-stained and Diff-Quik-stained smears can be used
- Gene mutations could be reliably detected down to 5% allele frequency
- Additional studies are underway to provide further analysis of cytology smear adequacy for ThyroSeq testing

## REFERENCES

1. Malapelle U, et al. Consistency and reproducibility of next-generation sequencing and other multigene mutational assays: A worldwide ring trial study on quantitative cytological molecular reference specimens. *Cancer Cytopathol.* 2017 Aug;125(8):615-626.
2. Nikiforova MN, et al. Analytical Performance of the ThyroSeq v3 Genomic Classifier for Cancer Diagnosis in Thyroid Nodules. *Cancer.* 2018 Apr;124(8):1682-1690.